

Advanced Test Equipment Rentals - www.atecorp.com 800-404-ATEC (2832)

VAISALA www.vaisala.com

GM70 Hand-Held Carbon Dioxide Meter for Spot-Checking Applications



The Vaisala CARBOCAP® Hand-Held Carbon Dioxide Meter GM70 is the demanding professional's choice for hand-held carbon dioxide measurement. The meter consists of the indicator (center) and probe, used either with the handle (left) or pump (right).

Features/Benefits

- Proven Vaisala CARBOCAP® reliability
- Two optional sampling methods: diffusion or pump aspiration
- User-friendly meter with multilingual user interface
- Numerical and graphical display of measurements
- Data can be logged and transferred to PC via MI70 Link software
- Wide selection of measurement ranges
- Easy recalibration using the interchangeable probes
- Suitable for field checking of fixed CO₂ instruments
- Short warm-up time
- Compact and versatile

The Vaisala CARBOCAP® Hand-Held Carbon Dioxide Meter GM70 is a user-friendly meter for demanding spot measurements in laboratories, greenhouses and mushroom farms. The meter can also be used in HVAC and industrial applications, and as a tool for checking fixed CO₂ instruments.

The GM70 has a short warm-up time and is ready for use almost immediately. It has a menu-based interface, a graphical LCD display and data logging capability. The optional MI70 Link Windows® software in combination with a USB connection cable provides an easy way to handle data in a PC environment.

Vaisala CARBOCAP® technology

The GM70 incorporates the advanced Vaisala CARBOCAP® sensor that has unique reference measurement capabilities. The measurement

accuracy is not affected by dust, water vapor or most chemicals. The GM70 has a two-year recommended calibration interval.

Two sampling methods

The handle is for hand-held diffusion sampling. The GM70 pump enables pump-aspirated sampling from locations difficult to access otherwise. It is also ideal for comparisons with fixed CO_2 transmitters.

Interchangeable probes

The GM70 uses the same probes as Vaisala CARBOCAP® Carbon Dioxide Transmitter Series GMT220 and Modules Series GMM220. By plugging different probes into the handle or pump, the user can easily change the measurement range of the GM70.

The meter can also be used as a calibration check instrument for Vaisala's GM20 and series, GMT/ M220 fixed CO_2 instruments. GMP220 probes can even be adjusted by using the GM70 meter.

The GM70 has two probe inputs. Vaisala's relative humidity and dewpoint probes can also be used simultaneously with ${\rm CO_2}$ measurement.

Technical data

CO2 volume concentration measurement

Measurement ranges		
High concentrations	02 %	
short probe (GMP221)	0 3 %	
	0 5 %, 0 10 %, 0 20 %	
Low concentrations	0 2000 ppm	
long probe (GMP222)	0 3000 ppm, 0 5000 ppm,	
	0 7000 ppm, 0 10,000 ppm	
Accuracy (including repeatability	non-linearity and calibration	
uncertainty) at 25 °C and 1013 hP	'a	
GMP221	$\pm (1.5\% \text{ of range} + 2\% \text{ of reading})$	
(applies for concentrations above 2% of full scale)		
GMP222	$\pm (1.5\% \text{ of range} + 2\% \text{ of reading})$	
Temperature dependence, typical	-0.3 % of reading / °C	
Pressure dependence, typical	+0.15% of reading/hPa	
Long-term stability	<±5 %FS/2 years	
Response time (63 %)		
GMP221	20 seconds	
GMP222	30 seconds	
Warm-up time	30 seconds, 15 minute full	
	specifications	

Measurement environment

Temperature	-20+60 °C (-4+140 °F)
Relative humidity	0 100 %RH non-condensing
Operation pressure	700 1300 hPa
Flow range (diffusion sampling)	0 10 m/s

Probe, handle & pump general

Housing material GMP221/222 probe GMH70 handle GM70 Pump Storage temperature Storage humidity GMH70 with GMP221/222 probe GM70 Pump with GMP221/222	Sensor	Vaisala CARBOCAP®	
GMH70 handle ABS/PC blend GM70 Pump aluminium casing Storage temperature -30+70 °C (-22+158 °F) Storage humidity 0100 %RH non-condensing Weight GMH70 with GMP221/222 probe GM70 Pump with GMP221/222	Housing material		
GM70 Pump aluminium casing Storage temperature -30+70 °C (-22+158 °F) Storage humidity 0100 %RH non-condensing Weight GMH70 with GMP221/222 probe GM70 Pump with GMP221/222	GMP221/222 probe	PC plastic	
Storage temperature Storage humidity Storage humidity Weight GMH70 with GMP221/222 probe GM70 Pump with GMP221/222 GM70 Pump with GMP221/222	GMH70 handle	ABS/PC blend	
Storage humidity 0 100 %RH non-condensing Weight GMH70 with GMP221/222 probe GM70 Pump with GMP221/222	GM70 Pump	aluminium casing	
Weight GMH70 with GMP221/222 probe 230 g GM70 Pump with GMP221/222	Storage temperature	-30+70 °C (-22+158 °F)	
GMH70 with GMP221/222 probe 230 g GM70 Pump with GMP221/222	Storage humidity	0 100 %RH non-condensing	
GM70 Pump with GMP221/222	Weight		
•	GMH70 with GMP221/222 probe	230 g	
probe 700 g	GM70 Pump with GMP221/222		
	probe	700 g	

MI70 Indicator general

MI/O indicator genera	ı
Menu languages	English, Chinese, French, Spanish,
	German, Japanese, Russian,
	Swedish, Finnish
Display	LCD with backlight, graphic trend
	display of any parameter,
	character height up to 16 mm
Max. no. of probes	2
Power supply	Rechargeable NiMH battery pack
	with AC-adapter
Analog output	0 1 VDC
Output resolution	0.6 mV
PC interface	MI70 Link software with USB or
	serial port cable
Data logging capacity	2700 points
Alarm	Audible alarm function
Operating temperature range	-10+40 °C (+14+104 °F)
Operating humidity range	non-condensing
Housing material	ABS/PC blend
Housing classification	IP54
Weight	400 g

Battery operation time

Continuous use	
with handle	better than 8h at $+20$ °C ($+68$ °F)
with pump	better than 5h at +20 °C (+68 °F) without load
Data logging use	up to a month, depending on logging interval

Electromagnetic compatibility

EN 61326-1, Portable Equipment.

Accessories

Accessories		
Connection cable for fixed CO ₂ instruments	_	
GMT220, GMM220, GMD20 and		
GMW20	GMA70	
MI70 Link software with USB cable	219687	
MI70 Link software with serial port	MI70LINK	
cable		
Analog output cable for 0 1 VDC	27168ZZ	
Calibration adapter	26150GM	
Carrying case	MI70CASE	
Battery, NiMH 4.8V	26755	
Spare probe	GMP221, GMP222	
(use the order form to define measurement range etc.)		
Nafion Membrane Tubing	212807GM	



For more information, visit us at sales@vaisala.com